

Molecular Weight Determinations of Heparin and Low-Molecular-Weight Heparins Using BioCore SEC-HP Columns

Application #: EN004-1

Applicable field: Heparin analysis

Author: Heparin analysis team at Asnail

Last modified date: May 31, 2019

Abstract

In this article, we describe the method for molecular weight determinations of unfractionated heparin and low molecular weight heparins (LMWHs) using the tailored BioCore SEC-HP columns. The results obtained show better resolution and longer durability than SEC columns widely used in the last decade.

Introduction

Molecular weight (MW) and its distribution are important characteristics of heparin and low molecular weight heparins (LMWHs), a family of anticoagulants with polydisperse, heterogeneous polysaccharide chains.

Current compendial methods for heparin molecular weight determinations involve separation on size-exclusion chromatography (SEC) columns, followed by detection with refractive index detectors (RID). While some SEC columns have been widely used in heparin industry before, they are subject to short life time.

BioCore SEC-HP columns, a series of SEC columns tailored for heparin MW determinations, are used for MW determinations of heparin and LMWHs. The results show good performance which meets the system suitability required by the latest methods specified in pharmacopoeias. Besides, these columns show very long life time, without performance loss even after 100 injections, which can reduce cost per injection for heparin industry.

Experiment

Equipment

HPLC system: Shimadzu LC-2030C

Detector: RID-20A refractive index detector

Columns

Table 1. SEC columns

| Cat. # | Description | Usage |
|-----------------------|--|-------------------------------|
| B213-050030-07830S-HP | BioCore SEC-300HP Analytical Column, 300 Å, 5 µm, 7.8×300 mm | Unfractionated heparin |
| B213-050050-07830S-HP | BioCore SEC-500HP Analytical Column, 500 Å, 5 µm, 7.8×300 mm | |
| B213-050015-07830S-HP | BioCore SEC-150HP Analytical Column, 150 Å, 5 µm, 7.8×300 mm | Low-molecular-weight heparins |
| B213-050015-04601S-HP | BioCore SEC-150 Guard Column, 5 µm, 4.6×10 mm | |

Reagents

Reference standard: Heparin low-molecular-mass for calibration CRS (European Directorate for the Quality of Medicines, Cat. # H019000), Heparin molecular weight calibrant (National Institutes for Food and Drug Control, Cat. # 140819), Heparin reference standard (National Institutes for Food and Drug Control, Cat. # 140818)

Sample: Low-molecular-weight heparin (provided by a customer)

Chromatographic conditions

Table 2. Chromatographic conditions

| Parameter | Value |
|-----------------------------|--|
| Mobile phase | 100 mM ammonium acetate |
| Elution mode | Isocratic |
| Flowrate | 0.5 mL/min |
| Injection volume | 25 µL (5 µL) |
| Column temperature | 30 °C (for system suitability test of UFH) 35 °C (for system suitability test of LWMHs and durability test) |
| Detector temperature | 30 °C (for system suitability test of UFH) 35 °C (for system suitability test of LWMHs and durability test) |
| Sample concentration | 5 mg/mL in mobile phase (for system suitability test of UFH) 10 mg/mL in mobile phase (for system suitability test of LMWHs) 10 mg/mL in water (for durability test) |

Procedure

The system suitability test of unfractionated heparin follows the method specified in Chinese Pharmacopoeia^[1].

The system suitability test of low-molecular-weight heparins and the durability test follow the method specified in European Pharmacopoeia^[2].

Results and discussion

System suitability test on heparin

Inject the heparin molecular weight calibrant solution into the HPLC system, and record the chromatogram, as shown in Figure 1.




Figure 1. Chromatogram of heparin molecular weight calibrant, analytical columns: BioCore SEC-300HP + SEC-500HP, guard column: none, injection volume: 25 μ L

By looking up the broad table of heparin molecular weight calibrant, we can identify each component peak corresponding to a specific molecular weight, as in Table 3.

Table 3. Peaks of heparin molecular weight calibrant

| Peak # | RT (min) | Response (mV) |
|--------|----------|---------------|
| 1 | 21.86667 | 40000 |
| 2 | 22.53333 | 36000 |
| 3 | 23.21667 | 32000 |
| 4 | 24.00000 | 28000 |
| 5 | 24.45833 | 26000 |
| 6 | 24.91667 | 24000 |
| 7 | 25.41667 | 22000 |
| 8 | 25.95833 | 20000 |
| 9 | 26.55000 | 18000 |
| 10 | 27.20000 | 16000 |
| 11 | 27.93333 | 14000 |
| 12 | 28.75833 | 12000 |
| 13 | 29.71667 | 10000 |
| 14 | 30.84167 | 8000 |
| 15 | 32.25000 | 6000 |




Figure 2. Calibration curve of heparin molecular weight calibrant

The calibration curve relating molecular weight to retention time can be obtained using a 3rd order polynomial fit, as shown in Figure 2.

The calibration formula is $y = -2.96724 \times 10^{-5}x^3 + 0.00147333x^2 - 0.093173x + 6.24575$, with a correlation coefficient $R = 0.999997$.

Inject the heparin reference standard solution into the HPLC system, and record the chromatogram, as shown in Figure 3.




Figure 3. Chromatogram of heparin reference standard, analytical columns: BioCore SEC-300HP + SEC-500HP, guard column: none, injection volume: 25 μ L

Table 4. Peaks of heparin reference standard

| Peak # | RT (min) | Molecular weight | Response (mV) |
|--------|----------|------------------|---------------|
| 1 | 24.22500 | 40000 | 1.177 |
| 2 | 24.46667 | 36000 | 1.383 |
| 3 | 24.75833 | 32000 | 1.712 |
| 4 | 25.14167 | 28000 | 2.225 |
| 5 | 25.39167 | 26000 | 2.608 |
| 6 | 25.66667 | 24000 | 3.071 |
| 7 | 26.00833 | 22000 | 3.663 |
| 8 | 26.39167 | 20000 | 4.339 |
| 9 | 26.83333 | 18000 | 4.974 |
| 10 | 27.35000 | 16000 | 5.312 |
| 11 | 27.93333 | 14000 | 5.103 |
| 12 | 28.59167 | 12000 | 4.379 |
| 13 | 29.35833 | 10000 | 3.333 |
| 14 | 30.26667 | 8000 | 2.127 |
| 15 | 31.43333 | 6000 | 0.992 |

The weight-averaged molecular weight (\bar{M}_w) can be calculated using the following formula

$$\bar{M}_w = \frac{\sum RI_i \times M_i}{\sum RI_i} = 16,611$$

Where RI_i is the RI response of the i^{th} peak, and M_i is the corresponding molecular weight calculated using the calibration formula.

The labelled weight-average molecular weight of heparin reference standard is 16,200. The difference between the labelled value and the calculated value is 411, which is within the range required by Chinese Pharmacopoeia, i.e. ± 500 .

To conclude, the correlation coefficient and the \bar{M}_w meet the system suitability requirements specified by the method of molecular weight determination in Chinese Pharmacopoeia.

System suitability test on low-molecular-weight heparins

SEC chromatograms of low-molecular-weight heparins, such as enoxaparin sodium, show fine structure of sub-peaks. The resolution between sub-peaks is an indicator of column performance.

Inject the heparin low-molecular-mass for calibration CRS solution into the HPLC system. The chromatogram of calibration CRS solution is shown in Figure 4.




Figure 4. Chromatogram of heparin low-molecular-mass for calibration CRS, analytical columns: BioCore SEC-150HP, guard column: BioCore SEC-150, injection volume: 25 μ L

Good resolution between sub-peaks of lower MW components can be observed, as shown in Figure 5. In particular, baseline separation is almost achieved between the two sub-peaks corresponding to MW 1,200 and 600.




Figure 5. Enlarged clip of Figure 4, showing sub-peaks

By looking up the broad standard table in the leaflet of heparin low-molecular-mass for calibration CRS, we can identify all peaks, as shown in Table 5.

Table 5. Peak identification corresponding to Figure 4

| Peak # | RT (Min) | MW |
|--------|----------|-------|
| 1 | 12.65000 | 18000 |
| 2 | 12.87500 | 15600 |
| 3 | 13.14167 | 13600 |
| 4 | 13.40833 | 12000 |
| 5 | 13.93333 | 9600 |
| 6 | 14.29167 | 8400 |
| 7 | 14.70000 | 7200 |
| 8 | 14.93333 | 6600 |
| 9 | 15.17500 | 6000 |
| 10 | 15.44167 | 5400 |
| 11 | 15.71667 | 4800 |
| 12 | 16.01667 | 4200 |
| 13 | 16.33333 | 3600 |
| 14 | 16.66667 | 3000 |
| 15 | 17.05000 | 2400 |
| 16 | 17.48333 | 1800 |
| 17 | 18.06667 | 1200 |
| 18 | 19.05833 | 600 |




Figure 6. Calibration of heparin low-molecular-mass for calibration CRS

The calibration curve relating molecular weight to retention time can be obtained using a 3rd order polynomial fit, as shown in Figure 6.

The calibration formula is $y = -0.00306062x^3 + 0.131503x^2 - 2.06044x + 15.4539$, with a correlation coefficient $R = 0.999579$.

Column durability

Heparin samples are very complex and cause great challenge for column durability. To test the durability of BioCore SEC-HP columns, we overlay chromatograms after multiple injections of the sample solution.

As we can see in Figure 7, BioCore SEC-HP column brings minimal variation between 30 injections of heparin samples, even without using the guard column.




Figure 7. Overlay of first 16 consecutive chromatograms of 30 injections of sample solution, injection volume 5 µL, obtained with BioCore SEC-150HP, no guard column

By incorporating the guard column, chromatograms obtained show almost no variation between 10 injections, even with the injection volume increased from 5 µL to 25 µL, as shown in Figure 8 and Figure 9.




Figure 8. Overlay of consecutive chromatograms of 10 injections of sample solution, injection volume 5 µL, obtained with BioCore SEC-150HP with BioCore SEC-150 guard column




Figure 9. Overlay of consecutive chromatograms of 10 injections of sample solution, injection volume 25 µL, obtained with BioCore SEC-150HP with BioCore SEC-150 guard column

Conclusion

As we can see in the last section, BioCore SEC-HP columns provide better resolution and longer durability for heparin MW determinations.

The superior performance of BioCore SEC-HP columns come in a large part from their monodisperse, low-absorption packing material.

To inspect the difference of packing materials, we take SEM of them.

As shown in Figure 10, the particle size distribution of BioCore SEC-HP is very narrow, with a mean size of 5 µm and CV<5%.




Figure 10. SEMs of particle size distributions of two packing materials

According to feedbacks from local heparin manufacturers, chromatographic profiles differ a lot after 10-20 injections if using competitors' columns. While BioCore SEC-HP columns demonstrate longer life time, with negligible performance lose after 50-60 injections of heparin samples.

The long life time of BioCore SEC-HP columns can be attributed to the superb mechanical strength of their packing material. To confirm that, we pack two columns under 6,000 psi using packing materials from BioCore SEC-HP columns and the other vendor's columns, followed by SEM of the packing materials. As shown in Fig. 9, no particle breakage occurs for BioCore SEC-HP packing material (a), while there is obvious breakage for the other (b).




Figure 11. SEMs of packing materials after 6,000 psi packing

A conclusion can be drawn that BioCore SEC-HP columns are very stable under high pressure and have great reproducibility in analyzing samples as complex as heparins. This durability performance gives a longer life time and can reduce cost.

Reference

1. Chinese Pharmacopoeia Commission. Chinese Pharmacopoeia (2015 Edition). Beijing: China Medical Science and Technology Press; 2015.6. p. Heparin Sodium; pp.517-518.
2. Council of Europe. European Pharmacopoeia. In: Supplement 9.8. Strasbourg: Council of Europe; 2019. p. Heparins, Low-Molecular-Mass; pp.6923-6924.